voco**Cimara****GB**

Instructions for use

CE 0482**Product description:**

Cimara is a ceramic repair material for the restoration of ceramic work. It contains the nanohybrid composite **GrandioSO**. Repairs of defective veneers can be carried out directly in the patient's mouth. The removal of dentures is unnecessary, as is the taking of impressions or the provision of temporary restorations. In many cases **Cimara** can also be used in situations where the metal frame is partially uncovered. The set contains all the materials required to create a durable, gap-free bond between ceramic and the composite.

Indications:

Intraoral repairs of defects in ceramic or composite veneer, fixed prosthetic work with metal frames.

Application:

Ensure that the working area is absolutely dry. The use of a rubber dam is recommended. A firm bond has to be created between composite and ceramic or metal frame in order to repair a defective ceramic veneer. Both surfaces must be prepared differently owing to their different material properties.

1. Preparing the metal surface:

Clean and roughen the alloy surface mechanically with a corundum grinding instrument. The adhesive strength of the bonding agent to be applied depends very heavily on the micro retentions created. If possible, make undercutts as well.

2. Preparing the ceramic margins:

Smooth defective ceramic veneer margins with a high speed (watercooled) bur and bevel to approximately 2 mm to achieve a firmer bond with the repair composite and better colour adaptation.

Drying the working area:

Dry ceramic and metal surfaces thoroughly with oil-free air.

3. Conditioning the ceramic margin:

Condition the bevelled ceramic margins with the enclosed **SiC Grinding Burs** at 6,000 - 10,000 rpm, without pressure, in a dry environment to create a clear surface free from grinding residues and ceramic debris.

It is essential to carry out this workstep thoroughly to ensure a firm bond. The surface of the **SiC Grinding Bur** absorbs the ceramic particles. Each grinding bur must therefore be used for one repair only.

4. Cleaning the bonding surface (metal, ceramic):

Carefully clean the bonding surfaces under dry conditions with the enclosed short-bristled brushes. It is essential that the cleaning of the working area and the subsequent worksteps are carried out under dry conditions.

5. Silanisation:

Open a capsule of **Cimara** coupling silane and dispense a few drops into a mixing dish. Apply a thin layer of **Cimara** coupling silane with a new brush onto the prepared ceramic surface and the exposed alloy. Allow the film to dry in the air for 2 minutes (do not use an air syringe) before carrying out the next step. The opened capsule must be used for one repair only and must not be kept for subsequent restorations (the coupling silane will decompose). Use only **Cimara** coupling silane to create an optimum bond.

6. Applying the Cimara opaquer on exposed metal frames:

Screw the enclosed metal cannula onto the **Opaquer LC** syringe. Apply a little **Opaquer LC** on the prepared metal surfaces and spread to a thin layer with a brush.

Then light-cure for 40 seconds. Conventional polymerisation devices (halogen or LED) are suitable for light-curing the material.

The light output should be at least 500 mW/cm². If the surface is not sufficiently covered, apply a second layer of **Opaquer LC** and light-cure again for 40 seconds.

7. Applying Cimara Adhesive:

Apply a thin layer of the **Cimara Adhesive** onto the prepared ceramic surface, distribute with a weak stream of air and light-cure for 20 seconds. With composite veneers, allow the material to take effect for 20 - 30 seconds. The inhibition layer formed on the surface must not be removed or contaminated because it is essential for creating a bond with the enclosed composite **GrandioSO**.

8. Repairing ceramic restorations with the nanohybrid composite GrandioSO:

The light-curing composite **GrandioSO** is available in 17 shades. The six most common shades are included in the set. Select the required shade and apply **GrandioSO** layer by layer. Polymerise each layer dependent on shade for 20 or 40 seconds with a lamp with a light output of at least 500 mW/cm². Then finish and polish with rotary instruments of the kind used in composite restoration techniques.

Note: Any areas of dental hard tissue involved in the repair must be treated with a dentine/dental enamel adhesive.

The Cimara system is also suitable for repairing composite veneers (omit workstep 3 for these type of repairs).

Information, precautionary measures:

To obtain the best and most durable results, ensure that no areas or only small areas of the metal frame are exposed. With larger metal surfaces (more than 50 %) and in stress zones, additional anchorage points should be created wherever possible (e. g. in the metal frame). If there is only a small area of ceramic surface available, or if there is no ceramic surface, it is recommended that the metal should be specially conditioned by sandblasting or that a secondary casting be fabricated.

Cimara contains methacrylates (Bis-GMA, UDMA, TEGDMA, hydroxyethyl methacrylate), amines, BHT and silane. **Cimara** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.

Reseal **Cimara Adhesive** bottle and **Cimara Opaquer LC** syringe after use.

Storage:

Close package immediately after use and avoid exposure to direct light. Store at temperatures of 4 °C to 23 °C. Storing the coupling silane at 8 °C extends its shelf life. Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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VC 60 AA1196 E1 0916 99 © by VOCO